



FlowWeave[®] & FlowWeave Bioseal[®]

Woven Polyester Vascular Grafts

- Most modern weaving technique guarantees high burst strength and low dilatation
- Specially texturised yarn for superior tissue incorporation and healing characteristics
- Aldehyde- and isocyanate free Bioseal[®] impregnation using dehydrothermal crosslinked collagen for blood tight grafts
- Concentric crimping and a guideline for safe and exact positioning
- Extraordinary soft and pliable handling

FlowWeave® & FlowWeave Bioseal®

Woven Polyester Vascular Grafts

Ordering information

Quick and reliable delivery in Europe within 24 hours.

Call our central office: +49 (0) 74 71 / 922 - 0

FlowWeave Bioseal®

Length (cm)	Diameter (mm)	Catalogue No.
15	12	45ST1512
15	14	45ST1514
15	16	45ST1516
15	18	45ST1518
15	20	45ST1520
15	22	45ST1522
15	24	45ST1524
15	26	45ST1526
15	28	45ST1528
15	30	45ST1530
15	32	45ST1532
15	34	45ST1534
15	38	45ST1538
30	10	45ST3010
30	12	45ST3012
30	14	45ST3014
30	16	45ST3016
30	18	45ST3018
30	20	45ST3020
30	22	45ST3022
30	24	45ST3024
30	26	45ST3026
30	28	45ST3028
30	30	45ST3030
30	32	45ST3032
30	34	45ST3034
30	38	45ST3038
60	10	45ST6010
60	12	45ST6012
60	14	45ST6014
60	16	45ST6016
60	18	45ST6018
60	20	45ST6020
60	22	45ST6022
60	24	45ST6024
60	26	45ST6026
60	28	45ST6028
60	30	45ST6030

These items can also be provided uncoated on request.

JOTEC®

SOLUTIONS FOR VASCULAR DISEASE

JOTEC GmbH

Lotzenäcker 23
D-72379 Hechingen
Tel. +49 (0) 74 71 / 922 - 0
Fax +49 (0) 74 71 / 922 - 100
www.jotec.net

JOTEC AG

Luzernerstr. 91
CH-5630 Muri
Tel. +41 (0) 56 670 91 - 00
Fax +41 (0) 56 670 91 - 02

